Vasopressin

¹²⁵I RIA Kit

For the quantitative determination of
human arginine vasopressin (AVP, ADH)

Instruction Manual

Manuel d'Instructions
Testanleitung
Manuale di Istruzioni

Catalog No./REF./KAT.-NR./Numero di Catalogo: 23065

DiaSorin

Stillwater, Minnesota 55082-0285, U.S.A.

TABLE OF CONTENTS

English	Page 1
Français	15
Deutsch	31
Italiano.....	46

VASOPRESSIN RADIOIMMUNOASSAY

1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE

This kit is intended for the quantitative determination of human arginine vasopressin (AVP, ADH) in plasma or other samples by radioimmunoassay (RIA) following an extraction using ODS-silica.

2. SUMMARY AND EXPLANATION

Arginine vasopressin (AVP) is a cyclic nine amino acid peptide. The structure of AVP is very similar to that of oxytocin differing in only two amino acids. Due to its powerful antidiuretic actions, vasopressin is also known as the antidiuretic hormone (ADH). Basically, AVP acts on the collecting tubule of the kidney, causing an increase in permeability to water and urea but not to larger molecules. AVP also possesses neurotransmitter functions, as well as releasing factor and peripheral humoral functions.¹

For many years, it has been known that AVP is produced in two locations in the hypothalamus, the supraoptic and the paraventricular nuclei.

AVP is transported from the supraoptic nucleus, via axons, through the neurohypo-physeal tract, to the posterior pituitary, where the axons end as terminal buttons on capillaries. It is from this region that secretion of AVP into peripheral circulation occurs, possibly in association with protein carrying molecules known as neurophysin.²⁻⁴

It seems likely that AVP from the paraventricular nucleus is secreted into capillaries of the hypothalamic hypophyseal portal system and conducted to the anterior pituitary where it may stimulate adrenocorticotropin (ACTH) release.

It is known that a number of both osmotic and non-osmotic stimuli will cause AVP release. A negative feedback system, which includes the kidney, central nervous system and thirst mechanism, exists and allows for regulation of body fluid concentration. AVP release is also influenced by a number of other factors including: emotional stress, posture, blood volume, temperature, and many pharmacological agents. Alcohol appears to weakly inhibit AVP secretion.⁵⁻⁷

The AVP assay is utilized in studies involving diabetes insipidus, syndrome of inappropriate ADH secretion (SIADH), syndrome of ectopic AVP production, and psychogenic water intoxication.^{5,8}

AVP RIA is very useful in separating central diabetes insipidus from nephrogenic diabetes insipidus. Central diabetes insipidus, which is marked by excessive thirst (polydipsia) and the passage of large amounts of dilute urine (polyuria), is caused by inadequate AVP production. Central diabetes insipidus can be either familial or idiopathic and can affect both males and females at any age. In contrast, nephrogenic diabetes insipidus is caused by failure of the renal tubules to respond to vasopressin, rendering the tubules unable to reabsorb water.

Nephrogenic diabetes insipidus, which is hereditary, generally affects males, and normally

SIADH refers to a release of vasopressin that is inappropriate to a low serum osmolality.^{11,5} It can be caused by a number of conditions including: pulmonary disease, head trauma (injury or surgical), local infection, and cancer. Certain neoplasms (especially bronchogenic carcinoma) produce AVP or AVP-like peptides. This ectopic vasopressin production syndrome is generally characterized by a sustained high elevation of AVP levels.¹²

Psychogenic water intoxication is caused by compulsive and excessive water consumption. It is not the result of a physical disorder; rather, as its name indicates, psychogenic water intoxication is a behavioral disturbance.

3. PRINCIPLE OF THE ASSAY

The DiaSorin AVP RIA consists of 2 procedures. The first involves extraction of the samples using the octadecasilyl-silica (ODS-silica) columns that are included in the kit. The extracted samples are then assayed using a disequilibrium RIA procedure to maximize sensitivity. There is a pre-incubation of sample and first antibody for 18 - 24 hours at 2-8°C; tracer is then added followed by a second incubation at 2-8°C. Phase separation is completed in 15-25 minutes at 20-25°C with a pre-precipitated complex of second antibody, carrier and PEG added in a single pipetting step.

4. REAGENTS PROVIDED IN THE KIT

Vasopressin 1% BSA-borate Buffer	1 vial/35 mL
Vasopressin Calibrator	1 vial/2.0 mL
Vasopressin Antiserum	1 vial/7 mL
¹²⁵ I Vasopressin	1 vial/7 mL
Precipitating Complex (GAGP-PPT)	1 vial/ 35 mL
Vasopressin Control Plasma	2 vials/2.5 mL
Number of tests	65

STORAGE: Upon receipt, and prior to reconstitution, store all reagents at 2-8°C. After reconstitution, store all reagents at -15° or lower until the expiration date on the label. Reagents should not be used past the expiration date.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

4.1 Vasopressin 1% BSA-borate Buffer: lyophilized reagent

1% BSA-borate buffer contains thimerosal as a preservative. Reconstitute the vial with 35 mL of purified water and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.2 Arginine Vasopressin Calibrator: lyophilized reagent

Purified arginine vasopressin, at a nominal concentration of 80 pg/mL, is made in BSA-borate buffer containing thimerosal. Exact values are assigned with each lot.

Reconstitute the vial with 2.0 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using. In order to obtain the entire calibrator curve, make serial dilutions by adding 500 μ L of calibrator to 500 μ L of 1% BSA-borate to give calibrators of 40, 20, 10, and 5 pg/mL. (These values will vary according to the exact concentration of the calibrator lot.)

Add 500 μ L of 80 pg/mL calibrator to 500 μ L of 1% BSA-borate and mix to give 40 pg/mL.

Add 500 μ L of 40 pg/mL calibrator to 500 μ L of 1% BSA-borate and mix to give 20 pg/mL.

Add 500 μ L of 20 pg/mL calibrator to 500 μ L of 1% BSA-borate and mix to give 10 pg/mL.

Add 500 μ L of 10 pg/mL calibrator to 500 μ L of 1% BSA-borate and mix to give 5 pg/mL.

The DiaSorin vasopressin calibrator has been calibrated against the World Health Organization (WHO) preparation 77/501. Any comparison of sensitivity with other products or procedures should be done with this reference standard. The calibrator demonstrates commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

4.3 Vasopressin Antiserum: lyophilized reagent

Guinea pig anti-vasopressin serum is diluted in BSA-borate buffer containing thimerosal. Reconstitute the vial with 7 mL of purified water and allow it to stand for 15 - 20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.4 ¹²⁵I Vasopressin: lyophilized reagent

(Lys⁹) Vasopressin is labeled with iodine-125 and diluted in BSA-borate-EDTA buffer containing thimerosal. Reconstitute the vial with 7 mL of purified water and allow it to stand for 15 -20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.5 Precipitating Complex (GAGP-PPT): lyophilized reagent

Normal guinea pig serum, pre-precipitated with goat anti-guinea pig serum and polyethylene glycol (PEG), is diluted in BSA-borate buffer with preservatives added. Reconstitute the vial with 35 mL of purified water; mix thoroughly until the suspension appears homogeneous and then allow it to stand for a minimum of 30 minutes at room temperature with occasional mixing.

4.6 Vasopressin Control Plasma: lyophilized reagent

Human plasma is spiked, if necessary, with the appropriate amount of arginine vasopressin to obtain a concentration within a specified range. Sodium azide is added as a preservative. Reconstitute each vial with 2.5 mL of purified water and allow them to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly and treat the quality control plasma as an unknown sample.

4.7 ODS-Silica: lyophilized reagent

Octadecasilyl-silica is enclosed in plastic columns (Sep-Paks). The ODS-silica columns included in the kit must be "primed" prior to applying sample (see Procedure section). The columns may be stored at room temperature.

5. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE

Not for internal or external use in humans or animals

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV), or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING THIMEROSAL

Some reagents in this kit contain thimerosal which contains a mercury compound. Disposal of elemental mercury, inorganic mercury, mercury oxides, and mercury compounds should be done in strict compliance with all local, state, and federal regulations.

WARNING: This product contains a chemical known to the State of California to cause cancer.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 1 μCi (37 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer, are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENT

- 6.1 The presence of abnormal particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

7. COLLECTION AND PREPARATION OF SPECIMEN

Two milliliters of EDTA plasma are required for the extraction in this procedure.

Collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. EDTA (7.2 mg/5 mL blood) is used as an anticoagulant. Centrifuge promptly for 15 minutes using 760 x g.* Separate the plasma from the cells and place in storage tubes. Vasopressin has been demonstrated to be stable at room temperature for several hours. For prolonged storage, freeze at -15°C or lower. All plastics, glassware, or other materials coming into contact with the specimen should be entirely free of any contamination.

NOTE: It has been observed that lower recoveries may be obtained when using heparinized plasma or serum. In addition, heparinized plasma will occasionally plug the Sep-Pak column.

8. EQUIPMENT AND MATERIAL REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Centrifuge to accommodate 12 x 75 mm tubes, temperature controlled.
- 8.3 Gamma scintillation counter capable of counting iodine-125.
- 8.4 Vortex.
- 8.5 Pipetting devices:
 - a. Micropipettors calibrated to deliver 100 µL and 200 µL.
 - b. Repeating dispensers calibrated to deliver 100 µL and 500 µL.
- 8.6 12 mL disposable syringes.
- 8.7 Methanol, reagent grade or better.
- 8.8 4% acetic acid.
- 8.9 1 M hydrochloric acid.
- 8.10 Compressed air or nitrogen supply.
- 8.11 37°C water bath.
- 8.12 Disposable glass tubes, 16 x 100 mm.
- 8.13 Purified water

9. EXTRACTION OF PLASMA WITH ODS-SILICA COLUMNS

- 9.1 Attach each Sep-Pak to a syringe hub.
- 9.2 Add 5 mL of methanol to each syringe.
- 9.3 Press through the syringe using the syringe plunger.
- 9.4 Remove the Sep-Paks.*
- 9.5 Remove the syringe plunger from the syringes.
- 9.6 Replace the Sep-Paks on syringes.
- 9.7 Add 10 mL of purified water to each syringe.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 9.8 Repeat steps 9.3 through 9.6.
- 9.9 Repeat step 9.7, then steps 9.3 through 9.6.
- 9.10 Acidify 2 mL of EDTA plasma or control plasma with 200 μ L of 1 M HCl. The pH should be 3 - 4. Load the sample in the syringe and push slowly, over a period of one minute, through the column.
- 9.11 Repeat steps 9.4 through 9.6.
- 9.12 Rinse the column with 20 mL of 4% acetic acid.
- 9.13 Repeat steps 9.3 through 9.6.
- 9.14 Collect the eluate in 16 x 100 mm tubes by eluting the peptide with 3 mL of methanol. Push slowly so that the methanol remains in contact with the ODS-silica for a minimum of 3 minutes.
- 9.15 Repeat steps 9.4 through 9.6.
- 9.16 Put through another 1 mL of methanol to ensure complete elution.
- 9.17 Evaporate the methanol eluate to dryness in a 37°C water bath using an air jet (compressed air) or nitrogen. Specimens must be completely dry. Do not leave in the 37°C water bath longer than required to dry.
- 9.18 Reconstitute the sample with 0.5 mL of 1% BSA-borate buffer. (A sample that is expected to be elevated can be diluted with more 1% BSA-borate at this step). Vortex and place at 37°C for ten minutes to ensure reconstitution of the dried sample. Vortex after incubation.
- 9.19 Assay 200 μ L in duplicate in the radioimmunoassay.

10. ASSAY PROCEDURE

- 10.1 Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures greater than 20-25°C. Mix all reagents gently before using.
- 10.2 Set up labeled 12 x 75 mm disposable tubes in duplicate according to the Scheme of the Assay, on the back page.
- 10.3 Add reagents as follows:
 - a. **Total count tubes**
Set aside until step 5
 - b. **Nonspecific binding tubes (NSB)**
200 μ L of 1% BSA-borate
 - c. **Calibrator 0**
200 μ L of 1% BSA-borate
100 μ L of vasopressin antiserum

- * Remove the Sep-Pak from the end of the syringe before removing the plunger to avoid contaminating the syringe contents with the fluid in the Sep-Pak.

d. Calibrators

200 μ L of vasopressin calibrator
100 μ L of vasopressin antiserum

e. Extracted control and unknown samples

200 μ L of extract
100 μ L of vasopressin antiserum

- 10.4** Vortex the tubes gently without foaming and incubate for 18-24 hours at 2-8°C.
10.5 Add 100 μ L of 125 I vasopressin to all tubes.
10.6 Vortex the tubes gently without foaming and incubate for 18-24 hours at 2-8°C.
10.7 Vigorously mix the GAGP-PPT and add 500 μ L to all the tubes except the total count tubes.
10.8 Vortex the tubes gently without foaming and incubate for 15-25 minutes at 20-25°C.
10.9 Centrifuge the tubes for 20 minutes using 760 x g^{*} at 20-25°C.
10.10 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a minimum time of 2 minutes. Blot the tubes with absorbent paper to remove any drops of supernatant that may be remaining on the tubes before turning the tubes upright.
10.11 Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy. (See Limitations of the Procedure section.)

11. PROCEDURAL COMMENTS

- 11.1** Assay all samples in duplicate to ensure confidence in values obtained.
11.2 Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
11.3 Some manufacturers' disposable borosilicate glass tubes yield elevated nonspecific bindings.
11.4 If you choose to aspirate the supernatant from the precipitate, be careful not to disturb the precipitate.
11.5 To completely monitor the consistent performance of an RIA there are additional factors which may be checked. DiaSorin suggests a check of the following parameters to assure consistent kit performance.

a. Total Counts

b. Maximum Binding

Average counts per minute (CPM) of Calibrator 0 Tube / Average CPM of Total Count Tubes.

$$* g = (1118 \times 10^8) (\text{radius in cm}) (\text{rpm})^2$$

c. Nonspecific Binding

Average CPM of NSB Tube / Average CPM of Total Count Tubes.

d. Slope of Calibrator Curve

For example, monitor the 80, 50, and 20% suppression points of the calibrator line.

12. QUALITY CONTROL

Each laboratory should include at least two control samples in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten assays. An acceptable range of values may then be obtained for these controls using ± 2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the control included with this kit. The range is printed on the control vial.

13. CALCULATION OF RESULTS

There are many methods in existence for calculating results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be neither linear or logarithmic scale. Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin quality control laboratory is % B/B₀ versus log concentration.

13.1 Calculate the average CPM for each calibrator, control and unknown sample.

13.2 Subtract the average CPM of the NSB tubes from all counts.

13.3 Divide the corrected CPM of each calibrator, control or unknown sample by the corrected CPM of the Calibrator 0.

$$\frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of Calibrator 0} - \text{CPM of NSB}} \times 100$$

13.4 Using two cycle semi-log or log-logit graph paper, plot percent B/B₀ for the vasopressin calibrators (vertical axis) versus the concentration (horizontal axis).

13.5 Draw a best-fit line through the points.

13.6 Interpolate the levels of vasopressin in the unknown samples from the plot.

13.7 Calculate the vasopressin value appropriately. For example: if 2 mL of plasma was extracted and reconstituted to 0.5 mL, divide by four; if 2 mL of plasma was extracted and reconstituted to 1.0 mL, divide by two.

13.8 If any sample reads greater than the highest calibrator it should be re-extracted with a lesser volume of plasma and/or reconstituted with more volume of 1% BSA borate.

13.9 If any extract should remain, it could also be diluted further and reassayed.

13.10 Calculate maximum binding by dividing CPM of Calibrator 0 by the average total counts obtained in the total count tubes.

TABLE II
DiaSorin Vasopressin RIA Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent Bound (B/T)	Percent (B/B ₀)	Graph Conc. (pg/mL)	Final Conc. (pg/mL)
Total Count	5,841 5,992	5,917					
NSB	192 205	199		3.4			
Calibrator 0	2,970 2,774	2,872	2,673	48.5	100		
Calibrators (pg/mL)							
A(5)	2,537 2,595	2,566	2,367		88.6		
B(10)	2,374 2,274	2,324	2,125		79.5		
C(20)	1,772 1,785	1,779	1,580		59.1		
D(40)	1,200 1,268	1,234	1,035		38.7		
E(80)	819 818	819	620		23.2		
Unknown Samples							
1	2,591	2,556 2,521	2,357		88.2	5.5	1.4
2	1,265 1,329	1,297	1,098		41.1	36.5	9.1

Typical sample data and a calibrator curve are shown in TABLE II and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

VASOPRESSIN SAMPLE CALIBRATOR CURVE


FIGURE 1

Reduction Data

The DiaSorin QC lab uses a smoothed spline curve fit.

14. LIMITATIONS OF THE PROCEDURE

14.1 Care must be taken so that all samples are extracted uniformly.

14.2 Counting times should be sufficient to prevent the introduction of statistical error in counting (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).

15. EXPECTED VALUES

Normal Range

Each laboratory should establish its own normal range. In male and female specimens (n = 27), DiaSorin determined the mean vasopressin value to be 1.4 pg/mL, with a normal range of ND-4.7 pg/mL. Normal values were calculated by extracting 2 mL of EDTA plasma. The final answers were divided by four to give the normal range of ND-4.7 pg/mL. The extraction and concentration step allows reporting of values that would otherwise be less than the sensitivity of the assay.

The values obtained by DiaSorin are consistent with those found by other researchers. In a study using thirty-one normal subjects, the normal range of plasma vasopressin was determined to be 1.4 ± 1.0 pg/mL.¹

16. SPECIFIC PERFORMANCE CHARACTERISTICS

16.1 Precision

Intra Assay Variation (values = pg/mL)

Sample Number	Mean Value	S.D	%C.V.
Low	4.9	0.54	11.2
Medium	22.7	1.38	6.1
High	41.5	2.25	5.4

Inter Assay Variation (values = pg/mL)

Sample Number	Mean Value	S.D	%C.V.
Low	5.0	0.14	2.8
High	40	1.78	4.4

16.2 Trueness: The assay trueness has been verified by the dilution test and the recovery test.

Linearity (Parallelism)

SERIAL DILUTION STUDY OF THREE PATIENT SAMPLES (VALUES = PG/ML)

Sample Number	Undiluted	1/2	1/4
1	5.4	5.5	4.7
2	9.3	8.9	9.0
3	3.5	3.1	2.8

Recovery (Accuracy) (Values = pg/mL)

Background	Calibrator Added	Expected Value	Measured Value	Percent Recovery	
Set No. 1	2.2	3.0	5.2	4.8	92
	2.2	5.0	7.2	6.3	88
	2.2	10	12.2	11	90
Set No. 2	1.9	3.0	4.9	5.0	102
	1.9	10	11.9	10.1	85
	1.9	20	21.9	20.3	93

16.3 Analytical Sensitivity

When defined as the apparent concentration at 3 standard deviations from the counts at maximum binding, the minimum detectable amount is 2.5 pg/mL. Any patient sample that reads less than 2.5 pg/mL before correcting for the dilution factor, cannot be distinguished from zero.

16.4 Analytical Specificity

Comparison of the reactivity of the vasopressin antibody was made with the following peptides:

Peptide	% Cross-reactivity
Arginine vasopressin	100
Lysine vasopressin	600
Oxytocin	<0.01
Vasotocin	0.14

SEE LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow reagents to reach above 25°C. Identify tubes in duplicate.
2. Extract controls and unknowns according section nine.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	NSB	Cal 0	Cal 1-5	Controls and unknown samples
1% BSA-borate Buffer	-	200 µl	200 µL	-	-
Calibrators (1-5)	-	-	-	200 µL	-
Extracted Controls	-	-	-	-	200 µL
Extracted Unknown	-	-	-	-	200 µL
Samples	-	-	100 µL	100 µL	100 µL
Antiserum					

4. Cover the tubes with parafilm and vortex gently.
5. Incubate for 18–24 hours at 2-8°C.
6. Dispense 100 µL of the tracer into all wells.
7. Cover and vortex gently.
8. Incubate for 18-24 hours at 2-8°C.
9. Dispense 500 µL of the GAGP-PPT into all wells except the total count tubes.
10. Incubate for 15-25 minutes at 20-25°C.
11. Centrifuge using 760 x g* for 20 minutes.
12. Decant the supernatants.
13. Count each tube in a gamma counter for 60 seconds or longer.

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$$*g = (1118 \times 10^6) (\text{raggio in cm}) (\text{rpm})^2$$

REFERENCES/BIBLOGRAPHIE/LITERATUR/BIBLIOGRAFIA

1. Cowley, A.W. Jr., W.C. Cusham, E.W. Quillen, Jr., M.M. Skelton and H.G. Langford, "Vasopressin Elevation in Essential Hypertension and Increased Responsiveness to Sodium Intake," **Hypertension**, 3:93, (1981).
2. Clark, G., P. Wood, L. Merrick and D.W. Lincoln, "Opiate Inhibition of Peptide Release from the Neurohumoral Terminals of Hypothalamic Neurones," **Nature**, 282:746, (1979).
3. Malvin, R.L., "Possible Role of the Renin-angiotensin System in the Regulation of Antidiuretic Hormone Secretion," **Federation Proceedings (Am. Physiological Soc.)**, 30:1383, (1971).
4. Pullan, P.T., B.H. Clappison and C.I. Johnston, "Plasma Vasopressin and Human Neurophysins in Physiological and Pathological States Associated with Changes in Vasopressin Secretion," **Journal of Clinical Endocrinology and Metabolism**, 49:580, (1979).
5. Lester, M.C. and P.B. Nelson, "Neurological Aspects of Vasopressin Release and the Syndrome of Inappropriate Secretion of Antidiuretic Hormone," **Neurosurgery**, 8:735, (1981).
6. Schrier, R.W. and J.P. Goldberg, "The Physiology of Vasopressin Release and the Pathogenesis of Impaired Water Excretion in Adrenal, Thyroid, and Edematous Disorders," **The Yale Journal of Biology and Medicine**, 53:525, (1980).
7. Weitzman, R.E., L. Farnsworth, R. MacPhee, C.C. Wang and C.M. Bennett, "Effect of Opposing Osmolar and Volume Factors on Plasma Arginine Vasopressin in Man," **Mineral and Electrolyte Metabolism**, 1:43, (1978).
8. Shimizu, K. and M. Hoshino, "Application of Vasopressin Radioimmunoassay to Clinical Study: Role of Vasopressin Hyponatremia and Hypernatremia and Some Other Disorders of Water Metabolism," **Contr. Nephrology**, 9:42, (1978).
9. Baylis, P.H. and G.L. Robertson, "Case Reports: Vasopressin Function in Familial Cranial Diabetes Insipidus," **Postgraduate Medical Journal**, 57:36, (1981).
10. Block, L.H., J. Furrer, R.A. Locher, W. Siegenthaler and W. Vetter, "Changes in Tissue Sensitivity to Vasopressin in Hereditary Hypothalamic Diabetes Insipidus," **Klinische Wochenschrift**, 59:831, (1981).
11. Haas, M., and S.M. Glick, "Radioimmunoassayable Plasma Vasopressin Associated with Surgery," **Archives of Surgery**, 113:597, (1978).
12. Yamaji, T., M. Ishibashi and S. Katayama, "Nature of the Immunoreactive Neurophysins in Ectopic Vasopressin-producing Oat Cell Carcinomas of the Lung," **Journal of Clinical Investigation**, 68:388, (1981).
13. Preibisz, J.J., J.E. Sealy, J.H. Laragh, R.J. Cody and B.B. Weksler, "Plasma and Platelet Vasopressin in Essential Hypertension and Congestive Heart Failure," Supplement I, **Hypertension**, 1:129, (1983).

Symbols used with devices

IVD

LOT

Ab

Ab

PEG

Ag

¹²⁵I

BUF

CAL

CONTROL

English	Français	Deutsch	Italiano
European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformità europea
Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Data di scadenza
Manufacturer	Fabricant	Hersteller	Fabbricante
Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso
In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnostica in vitro.
Lot No.	No. de lot	Chargen-Nr.	Lotto n°
Temperature limitation.	Limitation de température.	Temperaturbereich	Limite della temperatura
Antiserum	Antisérum	Antiserum	Antisiero
Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reagente precipitante
Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Tracciatore: antigene etichettato con ¹²⁵ I
Buffer	Tampon	Puffer	Tampone
Calibrator	Étalon	Kalibrator	Calibratore
Control serum	Sérum de contrôle	Kontrollserum	Siero di controllo
Radioactive	Radioactif	Radioaktiv	Radioattivo
Harmful	Nocif	Gesundheitsschädlich	Nocivo

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DiaSorin S.A.
11, rue Georges Besse
92160 Antony, France

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285 U.S.A.

Phone: 651-439-9710
FAX: 651-779-7847

For Customer Service in the U.S. and Canada Call Toll Free: 800-328-1482
In the United Kingdom Call: 44 118 9364200 FAX: 44 118 9792061
10397
27334 9/03

PRINTED IN U.S.A.