# 1,25-Dihydroxyvitamin D

<sup>125</sup>I RIA Kit

For the quantitative determination of  
1,25-Dihydroxyvitamin D in serum or EDTA plasma

## **Instruction Manual**

Manuel d'Instructions

Testanleitung

Manual de Instrucciones

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Bruksanvisning

**REF: 65100E**

**DiaSorin**

Stillwater, Minnesota 55082-0285, U.S.A.

**TABLE OF CONTENTS**

English .....	Page 1
Français .....	17
Deutsch .....	35
Español .....	53
Italiano.....	71
Svenska .....	88

## 1,25-DIHYDROXYVITAMIN D <sup>125</sup>I RIA KIT

### 1. INTENDED USE

FOR IN VITRO DIAGNOSTIC USE.

The 1,25-Dihydroxyvitamin D <sup>125</sup>I RIA is a competitive equilibrium radioimmunoassay intended for the quantitative determination of 1,25-dihydroxyvitamin D (1,25-(OH)<sub>2</sub>-D) in human serum or EDTA plasma to be used to assess 1,25-(OH)<sub>2</sub>-D deficiency associated with renal disease. Assay results should be used in conjunction with other clinical and laboratory data to assist the clinician in making individual patient management decisions in adult populations.

### 2. SUMMARY AND EXPLANATION

Vitamin D is derived from 2 sources: exogenous (dietary) and endogenous (biosynthesis, regulated by exposure to ultraviolet light). The exogenous or nutritional source includes foods containing naturally low levels of vitamin D<sub>2</sub> (e.g. milk, butter, cereals supplemented with vitamin D<sub>2</sub>), nutritional supplements in the form of readily available over-the-counter vitamins, and therapeutic formulations of the D<sub>2</sub> vitamins.<sup>1</sup> Vitamin D is not inherently active as it enters the circulation, by either dietary or photochemical routes. Biological activity is gained following a complex series of metabolic steps.<sup>2</sup>

It is now known that the metabolic activation of vitamin D is an intricately controlled process, subject to extensive alteration by variables including dietary calcium and phosphorus, degree of vitamin D deficiency, genetic deficiencies, parathyroid hormone concentrations, exposure to ultraviolet light, and degree of renal function.<sup>3</sup> The biosynthesis of the dihydroxylated forms of vitamin D<sub>3</sub> begins with the action of solar ultraviolet light on 7-dehydrocholesterol to form vitamin D<sub>3</sub> in the skin. Once vitamin D<sub>3</sub> enters the circulation it is rapidly taken up by the liver where it is metabolized to 25-hydroxyvitamin D<sub>3</sub> (25-OH-D<sub>3</sub>). The liver will also hydroxylate dietary vitamin D<sub>2</sub> to 25-hydroxyvitamin D<sub>2</sub> (25-OH-D<sub>2</sub>).<sup>2,4</sup>

Following hepatic hydroxylation, 25-OH-D is transported in association with vitamin D binding protein to the kidney where further hydroxylation takes place. The addition of a hydroxyl at position 1 yields 1,25-dihydroxyvitamin D (1,25-(OH)<sub>2</sub>-D). 1,25-dihydroxyvitamin D is the most potent naturally occurring vitamin D metabolite discovered so far, and its production is tightly regulated through concentrations of serum calcium, phosphorus, and parathyroid hormone. During times of calcium stress, 1,25-(OH)<sub>2</sub>-D is the most important vitamin D metabolite produced by the kidney.<sup>5,6</sup>

This is due to its essential role in the efficient active absorption of calcium and phosphorus, as well as their normal metabolism. Consequently, the measurement of 1,25-(OH)<sub>2</sub>-D is rapidly becoming an efficient tool in the research of diseases and conditions that affect the normal metabolism of phosphorus and calcium.<sup>7,8,9</sup>

### 3. PRINCIPLE OF THE ASSAY

The DiaSorin 1,25-(OH)<sub>2</sub>-D assay consists of a two-step procedure. The assay involves a preliminary extraction and subsequent purification of vitamin D metabolites from serum or EDTA plasma using C<sub>18</sub>OH cartridges.<sup>10</sup> Following extraction, the treated sample is then assayed using a competitive RIA procedure. The RIA method is based on a polyclonal antibody that is specific for both 1,25-(OH)<sub>2</sub>D<sub>2</sub> and 1,25-(OH)<sub>2</sub>D<sub>3</sub>. The sample, antibody and tracer are incubated for 2 hours at 20-25°C. Phase separation is accomplished after a 20 minute incubation at 20-25°C with a second antibody precipitating complex. After centrifugation and decantation, the bound fraction remaining in the pellet is counted in a gamma counter. Values are calculated directly from a calibrator curve of known concentrations. The final concentration of the 1,25-(OH)<sub>2</sub>D in serum and EDTA plasma samples is expressed as pg/mL.

#### 4. REAGENTS PROVIDED IN THE KIT

1,25-(OH) <sub>2</sub> D NSB BUFFER	1 vial / 3 mL
1,25-(OH) <sub>2</sub> D CALBRATOR O	1 vial / 20 mL
1,25-(OH) <sub>2</sub> D <sub>3</sub> CALIBRATOR 1-5	5 vials / 3 mL
1,25-(OH) <sub>2</sub> D ANTISERUM	1 bottle / 35 mL
<sup>125</sup> I 1,25-(OH) <sub>2</sub> D <sub>3</sub>	1 vial / 10 mL
1,25-(OH) <sub>2</sub> D <sub>3</sub> PRETREATMENT SOLUTION	1 bottle / 50 mL
1,25-(OH) <sub>2</sub> D GAR PRECIPITATING COMPLEX	2 vials / 35 mL
1,25-(OH) <sub>2</sub> D CONTROL SERUM	2 vials / 3 mL
95% ETHANOL	1 vial / 7 mL
Number of tests	100

**STORAGE:** Upon receipt, the kit should be stored at 2-8°C. After opening, store each reagent at 2-8°C until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

When reconstituting the contents of the vials, mix gently to avoid foaming. Store all reconstituted reagents at -15°C or lower immediately following use. Reagents from different batches must not be mixed.

**4.11,25-(OH)<sub>2</sub>D NSB Buffer:** Ready to use reagent

Potassium phosphate gelatin buffer containing ProClin® 300 (< 0.2%).

**4.21,25-(OH)<sub>2</sub>D 0 Calibrator:** Ready to use reagent

Phosphate buffer containing bovine serum proteins and ProClin 300 (< 0.2%).

**4.31,25-(OH)<sub>2</sub>D<sub>3</sub> Calibrators (1-5):** Lyophilized reagent

Five lyophilized (1,25-(OH)<sub>2</sub>D<sub>3</sub>) calibrators at concentrations ranging from 5 - 200 pg/mL containing human serum and ProClin 300 (<0.2%). Reconstitute each calibrator with 3.0 mL of calibrator zero. Exact concentrations are stated on the vial labels. The kit calibrators are calibrated using UV quantitation. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as recommended.

**4.41,25-(OH)<sub>2</sub>D Antiserum:** Ready to use reagent

Rabbit anti-1,25-(OH)<sub>2</sub>D serum diluted in phosphate-gelatin buffer containing ProClin 300 (<0.2%).

**4.5<sup>125</sup>I 1,25-(OH)<sub>2</sub>D<sub>3</sub>** Ready to use reagent

Radio-iodinated 1,25-(OH)<sub>2</sub>D<sub>3</sub> analogue diluted in an ethylene glycol phosphate buffer.

**4.61,25-(OH)<sub>2</sub>D<sub>3</sub> Pretreatment Solution:** Ready to use reagent

Potassium phosphate buffer.

**4.71,25-(OH)<sub>2</sub>D Controls: Level 1 (Normal), Level 2 (Elevated):** Ready to use reagent

Processed human serum pool, containing ProClin 300 (<0.2%), spiked with the appropriate amounts of 1,25-(OH)<sub>2</sub>D to obtain control concentrations within the specified ranges. Control 1 is within the normal range and Control 2 is within an elevated range. Expected ranges, as determined by DiaSorin, are stated on the vial labels.

**4.8Goat Anti-Rabbit (GAR) Precipitating Complex:** Lyophilized reagent

Normal rabbit serum, pre-precipitated with goat anti-rabbit serum and polyethylene glycol (PEG), is diluted in BSA-borate buffer with 0.1% sodium azide and other preservatives added (lyophilized). Reconstitute the vial with 35 mL of distilled or deionized water; mix thoroughly until the suspension appears homogeneous and then allow it to stand for a minimum of 30 minutes at room temperature with occasional mixing.

**4.995% Ethanol:** Ready to use reagent  
95% ethanol and 5% water.

**NOTE: C<sub>18</sub>OH cartridges are also required for this procedure. These cartridges must be ordered separately under DiaSorin Cat. No. 65101E.**

**5. WARNINGS AND PRECAUTIONS**

FOR IN VITRO DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

**CAUTION:** This device is to be received, acquired, possessed, and used only by physicians, clinical laboratories, hospitals or research facilities. Testing is to be performed only by properly qualified and trained laboratory personnel.

**REAGENTS CONTAINING HUMAN SOURCE MATERIAL**

**Treat as potentially infectious.**

Each serum/plasma donor unit used in the preparation of this product has been tested by an U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus, hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4<sup>th</sup> ed., May 1999 or current edition.

**REAGENTS CONTAINING SODIUM AZIDE**

**CAUTION:** Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, U.S. 1976.

**European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)**

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

**REAGENTS CONTAINING ETHANOL**

**European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)**

R11 - Highly Flammable

S16 - Keep away from sources of ignition - No smoking

S43 - In case of fire, use dry chemical or carbon dioxide

## REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 5.5  $\mu\text{Ci}$  (204 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**WARNING:** This product contains a chemical known to the State of California to cause cancer.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

## 6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 The presence of abnormal particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

## 7. COLLECTION AND PREPARATION OF SERUM AND PLASMA

The 1,25-Dihydroxyvitamin D kit requires 500  $\mu\text{L}$  of either serum or EDTA plasma to perform the assay extraction; a volume of 1.5 mL will permit repeat analysis and provide adequate pipetting volume as well.

Either human serum or plasma may be used in this kit. The anticoagulant EDTA may be used with this assay. A fasting specimen is recommended, but not required. Blood should be collected aseptically by venipuncture in a 5 or 10 mL evacuated glass tube.

For serum, allow the blood to clot at room temperature (15-25°C). Centrifuge for 15 minutes using approximately  $760 \times g^*$  to obtain hemolysis free sera. No additives or preservatives are required to maintain integrity of the sample. All plastics, glassware or other material coming into contact with the specimen should be entirely free of any contamination.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

Store serum or plasma samples at -15°C or lower. Samples should not be repeatedly frozen and thawed. However, a study performed by DiaSorin showed no significant change in values following 3 freeze-thaw cycles of the samples. Samples have been stored at DiaSorin up to 6 months at -15°C or lower with no significant change in results.

## 8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1 C<sub>18</sub>OH Cartridges (24 cartridges) ordered separately under DiaSorin Cat. No. 65101E.
- 8.2 Disposable borosilicate glass tubes, 12 x 75 mm and 13 x 100 mm.
- 8.3 Test tube rack.
- 8.4 Parafilm or equivalent for covering test tubes.
- 8.5 Foam rack for decanting or equivalent.
- 8.6 Centrifuge capable of accommodating 12 x 75 mm tubes and attaining 1800 x g\*.
- 8.7 Gamma counter capable of counting <sup>125</sup>I.
- 8.8 Vortex mixer.
- 8.9 Pipetting devices:
  - a. Micropipettors calibrated to deliver 75 µL, 300 µL and 500 µL.
  - b. Repeating dispensers capable of delivering 50 µL, 300 µL, and 500 µL.
  - c. Volumetric pipettes for reconstituting calibrators, 3.0 mL.
- 8.10 Deionized or distilled water.
- 8.11 Organic solvents:

The following solvents will be needed for the preparation and extraction of samples.

**Use only HPLC-grade solvents.**

**Do not expose any of the organic solvents to acid washed glassware or any other acidic conditions.**

  - a. Acetonitrile.
  - b. Methanol.
  - c. Hexane. (IMPORTANT: n-hexane percentage must be greater than 95%)
  - d. Methylene chloride (must **not** contain alcohol as a preservative).
  - e. Isopropanol (2-propanol).
- 8.12 Recommended device for the C<sub>18</sub>OH cartridges:
  - a. Vac Elut sample processing station. Processes 24 columns per run. Available from Analytichem International or DiaSorin. To order from DiaSorin, use DiaSorin Cat. No. 11610.
- 8.13 Materials needed for drying samples:
  - a. Nitrogen gas.
  - b. Drying manifold with 37°C (±2°C) heating block or water bath: N-EVAP or MULTI-VAP available through Organomation Associates (Berlin, MA) or Reacti-Vap II and Reacti-Therm III available through Pierce Chemical Company (Rockford, IL).

\*  $g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$

## 9. PREPARATION OF ORGANIC SOLVENTS FOR THE COLUMN EXTRACTION

**NOTE:** The Organic solvent mixtures should be prepared before starting the Preliminary Extraction Procedure.

### PRECAUTIONS:

Prepare solvent mixtures as described in TABLE I; a 500 mL volume preparation is recommended by DiaSorin. Larger or smaller volumes can be prepared as long as the proportions remain consistent. However, volumes should be large enough to minimize measurement error. Measure each solvent independently for best results.

### PRECAUTIONS:

Use distilled or deionized water where applicable.

**TABLE I**  
Solvent Preparation

Solvent Mixture	Materials	1 L	500 mL
70:30	methanol water	700 mL 300 mL	350 mL 150 mL
90:10	hexane methylene chloride	900 mL 100 mL	450 mL 50 mL
99:1	hexane IPA	990 mL 10 mL	495 mL 5 mL
92:8	hexane IPA	920 mL 80 mL	460 mL 40 mL

## 10. PRELIMINARY EXTRACTION PROCEDURE

- 10.1** Reconstitute each lyophilized calibrator with 3.0 mL of calibrator zero. Mix the calibrators well and let stand for 15-20 minutes to ensure complete reconstitution. Thaw any frozen reagents completely. Allow all reagents to reach room temperature. Do not allow reagents to reach temperatures above 25°C. Mix all reagents well before use.
- 10.2** Pipette 500 µL of each calibrator (0-5), control and patient sample into labeled 12 x 75 mm borosilicate glass tubes.
- 10.3** Add 500 µL of acetonitrile to each sample; vortex intermittently, at least 3 times, during a time period of 10 minutes.
- 10.4** Centrifuge tubes at 760 x g\* for 10 minutes at 20-25°C.
- 10.5** Decant the sample supernatants into labeled 12 x 75 mm or 13 x 100 mm (if preferred) test tubes; discard pellets.
- 10.6** Add 500 µL of Pretreatment Solution to each tube and vortex.  
**NOTE:** This kit is designed to accommodate up to 48 extractions including the calibrators and controls plus 40 unknowns. This means that 24 samples can be extracted on the VacElut immediately and then the remaining set of 24 samples can be extracted afterwards. The column extractions should be done as soon as possible after the pretreatment addition. However, once column extracted, samples may be stored for 96 hours at -15°C or below.
- 10.7** The samples are now ready to be applied to the C<sub>18</sub>OH cartridges.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

## 11. COLUMN PROCEDURE

- 11.1 For an explanation of the VAC ELUT, refer to the VAC ELUT User's Guide. Assemble the station as instructed.
- 11.2 Label a 13 x 100 mm borosilicate glass tube for each calibrator (0-5), control, and unknown sample. Place these tubes in the VAC ELUT. Make sure the VAC ELUT lid is in the "WASTE" position and the tubes are positioned to collect the desired eluates. Place the C<sub>18</sub>OH cartridges on to the VAC ELUT cover.
- 11.3 Add the solvent mixtures to the specified cartridge as described in TABLE II.

### IMPORTANT:

#### Vacuum Usage:

Turn on the vacuum after each solvent addition and allow the solvent mixture to flow **completely** through the cartridge before proceeding to the next solvent. Turn the vacuum on and off at the recommended steps in TABLE II (if preferred, the vacuum may be turned off between each solvent addition). The vacuum should be set at 10 inches (254 mm of Hg) or lower. Elute each solvent mixture through the "WASTE" outlet and into appropriate waste containers until the final collection step.

#### Sample Application:

The samples can either be decanted directly onto the cartridge or applied with a pipette.

Do not let the cartridges air dry for more than 5 minutes between applications.

#### C<sub>18</sub>OH Cartridge Pretreatment:

The C<sub>18</sub>OH cartridges should be washed with 5 mL of 90:10 (hexane: methylene chloride), 5 mL of IPA and 5 mL of methanol prior to first time usage. The 90:10 (hexane: methylene chloride) solvent mixture can be prepared by measuring 900 mL of hexane and 100 mL of methylene chloride. Place the new, unused cartridges on the VAC ELUT apparatus, set in the "WASTE" position, and add 5 mL of 90:10 (hexane: methylene chloride) followed by 5 mL of IPA to each column and then 5 mL of methanol. Allow each solvent mixture to pass entirely through the cartridges before proceeding to the next mixture. After this initial preparation, it will not be necessary to repeat this step since the cartridges will be regenerated during the extraction procedure.

**NOTE:** If additional samples are to be extracted (more than 24), regenerate the cartridges using the same method outlined in TABLE II. The cartridges can be re-used up to 30 times.

**NOTE:** The C<sub>18</sub>OH Cartridge material is held in place using a porous fiber frit. Discard cartridges with dislodged or missing frits as the C<sub>18</sub>OH material may be lost.

**TABLE II**

Description and Analysis of Extraction	Extraction Steps	Eluant Handling
<p><b>Column Preparation/Regeneration</b> Removes: Interfering substances from previous usage (step 1)</p>	<p>1. Add 1 mL Methanol, Turn on Vacuum. 2. Turn off Vacuum.</p>	Discard "WASTE"
<p><b>Sample Application</b></p>	<p>3. Apply all Samples, Turn on Vacuum. from Preliminary Extraction Procedure (step 10).</p>	Discard "WASTE"
<p><b>Purification/Removal of Vitamin D Metabolites</b>  Removes: Interfering polar lipids, salts, and pigments (step 4)  25(OH)D (Step 5)  Remaining 25(OH)D and 24,25(OH)<sub>2</sub>D/25,26(OH)<sub>2</sub>D (step 6)</p>	<p>4. Add 5 mL 70:30 Methanol/Water (deionized or distilled). 5. Add 5 mL 90:10 Hexane/Methylene Chloride. 6. Add 5 mL 99:1 Hexane/Isopropanol. 7. Turn off the Vacuum.</p>	Discard "WASTE"
<p><b>Collection of 1,25-(OH)<sub>2</sub>D</b>  Elution of purified 1,25-(OH)<sub>2</sub>D (step 9)</p>	<p><b>IMPORTANT!</b> 8. Switch Vac Elut to "Collect:"  9. Add 3 mL 92:8 Hexane/Isopropanol, Turn on the Vacuum.  10. Turn off the Vacuum.</p>	Discard "WASTE"

**12. DRYING AND RECONSTITUTION OF ELUATES**

- 12.1 Place the sample tubes containing the eluates in a heating block or water bath at 37°C (±2°C) for drying.
- 12.2 Dry the eluates under a hood using 2-4 psi of nitrogen gas (drying time 20-30 minutes).
- 12.3 Remove the tubes immediately after the eluate has dried.
- 12.4 Reconstitute each of the dried extracts with 50 µL of 95% ethanol; vortex gently at a low to medium speed. Add 125 µL of tracer into the same tubes containing the 50 µL of 95% ethanol, vortex gently again at a low to medium speed. **NOTE:** The vortex steps are very important for ensuring that the sample has been properly reconstituted and for obtaining good precision. **CAUTION:** Keep the vortex at the lower portion of the tube to prevent loss of sample volume which will ensure sufficient pipetting volume for the assay. Label one extra tube for Total count and one tube for NSB. Add 50 µL of 95% ethanol and 125 µL of tracer into both tubes. These will be used to create the TC (Total Count) and NSB duplicate tubes for the assay set up.
- 12.5 Perform the assay **immediately** after reconstituting the samples. For best results when performing larger assays, reconstitute 12 - 15 samples at a time and then pipette into the assay before reconstituting the next 12 - 15 samples.

### 13. ASSAY PROCEDURE

- 13.1 Set up labeled 12 x 75 mm disposable borosilicate glass tubes in duplicate for each calibrator (0-5), control and sample according to scheme of the assay. Carefully add 75  $\mu$ L of the reconstituted calibrator, control, and sample extracts into the duplicate assay tubes.  
**CAUTION:** This step must be performed carefully since there is only 25  $\mu$ L excess in each tube.
- 13.2 Refer to "Drying and Reconstitution of Eluates" step number 4. Add 75  $\mu$ L from the TC (Total Count) tube into duplicate assay tubes. Repeat this step for the NSB duplicate assay tubes as well.
- 13.3 Add 300  $\mu$ L of NSB buffer into the NSB tubes.
- 13.4 Add 300  $\mu$ L of primary antibody into all tubes except the Total Counts and NSB tubes.
- 13.5 Mix well; incubate for 2 hours ( $\pm$ 15 minutes) at 20-25°C.  
**NOTE:** Reconstitute the GAR Precipitating Complex with 35 mL of distilled or deionized water, mix thoroughly until the suspension appears homogenous and then allow it to stand for a minimum of 30 minutes before use at room temperature with occasional mixing before and during use.
- 13.6 Add 500  $\mu$ L of well-mixed GAR Precipitating Complex into all tubes except the Total Count tubes. Incubate for 20 minutes ( $\pm$ 5 minutes) at 20-25°C.
- 13.7 Centrifuge all tubes for 20 minutes at 20-25°C at 1800 x g\*, except the Total Count tubes.
- 13.8 Decant the supernatants, except the Total Count tubes, using a foam rack tube holder or equivalent by inverting the rack into an appropriate waste container. Place the inverted rack onto absorbent paper for 2 - 3 minutes. Blot the tubes gently to ensure all liquid is removed.
- 13.9 Measure the radioactivity by counting all tubes for 1 minute on a gamma counter. Tubes should be counted for a minimum of 1 minute (see Limitations of Procedure section).

### 14. PROCEDURAL COMMENTS

- 14.1 Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
- 14.2 Correct solvent proportions are critical for good recoveries. Prepare solutions in volumes large enough to minimize measurement error.
- 14.3 If any sample reads greater than the highest calibrator, it should be re-assayed by diluting with the kit calibrator zero prior to extraction. Results should be multiplied by the appropriate dilution factor. For example, mix 500  $\mu$ L of the sample with 500  $\mu$ L of the calibrator zero, then assay 500  $\mu$ L of this sample mixture as per the product insert.
- 14.4 No assay drift was observed in the average time it takes to perform a 100 tube assay.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

**14.5** In order for a laboratory to completely monitor the consistent performance of an RIA assay, there are additional factors which may be checked. DiaSorin suggests a regular check of the following parameters to assure consistent kit performance.

**a. Total Counts**

**b. Maximum Binding**

Average counts per minute (CPM) of 0 calibrator tubes/Average CPM of total count tubes.

**c. Nonspecific Binding**

Average CPM of the NSB tubes/Average CPM of total count tubes.

**d. Slope of the Calibrator Curve**

The 50% suppression can be monitored.

**15. QUALITY CONTROL**

Each laboratory should include at least two controls (one normal range control and one elevated range) control in every assay to monitor kit performance. Commercially available controls or the 2 reference controls provided with the kit may be utilized.

The controls used should be treated as unknown specimens and assayed in duplicate. Quality control charts should be maintained by the laboratory to follow performance of the controls. Acceptable performance limits should be determined by each individual laboratory for each level of control using statistically based methods designed to detect both random and systematic errors. Control results must meet the laboratory's criteria for acceptability prior to reporting patient test results.<sup>12,13,14</sup>

**16. CALCULATION OF RESULTS**

There are many methods in existence for calculating results of RIA assays. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either a linear or logarithmic scale. Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin Quality Control Laboratory is a %B/B<sub>0</sub> versus log concentration calculation method based on a Spline Smoothed curve fit program.

**16.1 Calculation of Percent B/B<sub>0</sub>**

- a. Calculate the average CPM for each calibrator, control and unknown sample.
- b. Subtract the average CPM of the NSB tubes from all counts.
- c. Divide the average corrected CPM of each calibrator, control or unknown sample by the average corrected CPM of the 0 calibrator and multiply by 100.

$$\frac{\text{avg. CPM (Calibrator or Unknown Sample)} - \text{avg. CPM (NSB)}}{\text{avg. CPM (0 Calibrator)} - \text{avg. CPM (NSB)}} \times 100 = B/B_0 (\%)$$

**16.2 Use of Calibrator Curve Plot**

- a. Using 2 cycle log-logit or semi-log graph paper, plot percent B/B<sub>0</sub> (%) for the 1,25 (OH)2D calibrators on the ordinate (Y axis) versus the calibrator concentration on the abscissa (X axis).

NOTE: Automated data reduction programs may also be utilized in the analysis of data. DiaSorin utilizes Multi-Calc (Pharmacia) with a LIN-LOG smooth-SPLINE fit program. Other data reduction methods must be validated before incorporating for regular use.

- b. Draw a best-fit line through the points.
- c. Interpolate the levels of 1,25-(OH)<sub>2</sub>D in the samples from the plot of the calibrators.
- d. If any unknown sample has been diluted, correct for the appropriate dilution factor.
- e. The reportable range of the assay is 5.0 pg/mL to 200 pg/mL. Any value that reads below the lowest calibrator, 5.0 pg/mL, is an extrapolated value and may be reported as "less than 5 pg/mL".
- f. Calculate maximum binding by dividing CPM of 0 calibrator by the average total counts obtained in the total count tubes.

Typical sample data for the 1,25-(OH)<sub>2</sub>D RIA is shown in TABLE IV and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

**TABLE IV**  
DiaSorin 1,25-(OH)<sub>2</sub>D RIA Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent Bound (B/T)	Percent (B/B <sub>0</sub> )	Conc. (pg/mL)
Total Count	40,792 40,578	40,685				
NSB	1,438 1,439	1,439		3.5		
0 Calibrator	18,006 18,119	18,063	16,624	44.4	100.0	
Calibrators (pg/mL)						
1 (5.0)	16,710 16,770	16,740	15,302		92.0	
2 (15.0)	14,941 14,916	14,929	13,490		81.2	
3 (30.0)	13,087 13,124	13,106	11,667		70.2	
4 (75.0)	9,820 9,765	9,793	8,354		50.2	
5 (200)	6,412 6,464	6,438	5,000		30.0	
Controls						
Level 1: normal range	13,275 13,530	13,403	11,964		72.0	27.3
Level 2: high range:	8,206 8,124	8,165	6,727		40.5	119

DiaSorin 1,25-(OH)<sub>2</sub>D RIA Sample Calibrator Curve


FIGURE 1

#### DATA REDUCTION

The DiaSorin QC lab uses a smoothed spline curve fit.

#### 17. LIMITATIONS OF THE PROCEDURE

- 17.1 Counting times should be sufficient to prevent the introduction of error due to counter inefficiency (for example, accumulation of 2,000 CPM will yield 5% error, 10,000 CPM will yield 1% error).
- 17.2 Assay results should be used in conjunction with other clinical and laboratory data to assist the clinician in making individualized patient management decisions.
- 17.3 The performance characteristics of this assay have not been established in a pediatric population.

#### 18. EXPECTED VALUES

##### Reference Interval for Normal Donors

It is important for each laboratory to establish its own reference interval representative of its typical population. However, 1,25-(OH)<sub>2</sub>D values were collected from 123 apparently healthy volunteers from three Midwestern U.S.A cities in a clinical trial conducted at three sites during the summer and autumn months. These 123 healthy donors were composed of a racially diverse set of 37 males and 86 females, with an age range of 21-68 years. The mean 1,25-(OH)<sub>2</sub>D value for the entire sample (n=123) was 43.9 pg/mL. The 95% reference interval estimated by a nonparametric method (following NCCLS guideline C28-A2<sup>15</sup>) was 25.1-66.1 pg/mL.

##### End-stage renal disease patients

A clinical trial was conducted to evaluate the levels of 1,25-(OH)<sub>2</sub>D that would be found in an adult population diagnosed with end-stage renal disease. A total of 87 such volunteers (49 males, 38 females, age range of 19-84) from three Midwestern U.S.A.

cities were collected and assayed. The observed range of 1,25-(OH)<sub>2</sub>D values for this sample (n=87) was 1.6-17.3 pg/mL. The 95% upper reference limit estimated by the nonparametric procedure is 14.2 pg/mL.

## 19. SPECIFIC PERFORMANCE CHARACTERISTICS

### 19.1 Precision

Assay precision was evaluated by DiaSorin based on the principles of the NCCLS Guideline (EP5-A).<sup>16</sup> Three human serum-based control levels with 1,25-(OH)<sub>2</sub>-D concentrations distributed across the assay range were assayed over 25 assay days, spanning more than 60 operating days. Multiple technicians, as well as multiple lot numbers for all components were included. The combined results were evaluated by analysis of variance (ANOVA) and are summarized in the following table.

Sample	N	Mean pg/mL	Within-run		Between Day		Total	
			S.D.	%CV	S.D.	%CV	S.D.	%CV
LOW	25	25.8	1.76	6.8	3.8	14.6	4.0	15.3
MID	25	41.3	3.16	7.7	4.6	11.1	5.1	12.3
HIGH	25	105.2	11.86	11.3	11.8	11.2	14.4	13.7

### 19.2 Trueness: The assay trueness has been checked by the linearity test and the recovery test.

#### Linearity (Parallelism)

Linearity of dilution was evaluated by DiaSorin consistent with the recommendations found in NCCLS guideline, EP6-P.<sup>17</sup> Three serum patient sample pools were prepared by serial dilution with calibrator zero and aliquots individually frozen at -20°C. Each sample and dilution were assayed in multiple replicates over three different assay dates. Expected values were determined by the undiluted values of each sample pool multiplied by the dilution factor. The data is summarized below and the composite results plotted as a linear regression of Expected vs. Observed Value.

Sample #	Dilution	Expected Value (pg/mL) (N=10)	Mean Observed Value (pg/mL) (N=10)
1	NEAT	49.6	49.6
	1:2	24.8	25.3
	1:4	12.4	11.1
	1:8	6.2	4.4
2	NEAT	50.2	50.2
	1:2	25.1	24.2
	1:4	12.6	12.2
	1:8	6.3	4.6
3	NEAT	62.0	62.0
	1:2	31	31.9
	1:4	15.5	13.4
	1:8	7.8	6.2

### Recovery

The ability of the assay to quantitatively recover all of the analyte present in clinical samples was evaluated by the addition of different amounts of freshly prepared pure 1,25-(OH)<sub>2</sub>-D antigen to three patient sample pools. Three concentrations were chosen and assayed in duplicate, and percent recovery determined. The mean % recovery by this method was 101%.

Sample #	Init. Conc (pg/mL) 1.0 mL	Spike Amount (pg)	*Expected (pg/mL)	Observed (pg/mL)	% Recovery
1	30.8	50.0	77.7	75.6	97
		62.5	88.9	87.2	98
		75.0	99.8	103	103
2	42.8	50.0	89.2	89.5	100
		62.5	100.3	98.9	99
		75.0	111.0	118	106
3	40.3	50.0	86.8	83.9	97
		62.5	97.9	103	105
		75.0	108.8	118	108

\*Expected concentration calculation includes dilution factor introduced by spiking solution.

### 19.3 Analytical Sensitivity

The sensitivity of this assay, when defined as the lowest quantity differentiated from zero at 2 calibrator deviations below the means cpm of the calibrator zero (n=20), has been shown to be  $\leq 2.0$  pg/mL.

### 19.4 Analytical Specificity

Data on the cross-reactivity of the anti serum used in this kit is expressed as the ratio of 1,25-(OH)<sub>2</sub>D<sub>3</sub> concentration to the cross-reacting substance concentration at 50% inhibition of maximum binding.

Analyte	Conc. At 50% B/Bo	% Cross Reactivity
1,25 D <sub>3</sub>	136 pg/mL	100
1,25 D <sub>2</sub>	125 pg/mL	109
24,25 D <sub>3</sub>	176,000 pg/mL	<0.08
25,26 D <sub>3</sub>	148,000 pg/mL	<0.08
25 D <sub>3</sub>	>>>1 mg/mL	<0.01

### 19.5 Interfering Substances

An interference study was conducted to determine whether elevated levels of common endogenous substances could negatively impact assay results. The evaluation was consistent with NCCLS guidelines (EP7-P).<sup>18</sup> Three human serum based samples containing low, medium or high levels of 1,25-(OH)<sub>2</sub>-D<sub>3</sub> were tested with either additions of the test substance or as controls spiked with identical volumes of test substance vehicle. Each sample was extracted twice generating four replicates. Results presented below demonstrate that there was no interference by any of the tested substances, as determined by statistical testing by ANOVA (95% confidence interval), or by the mean value of the spiked samples exceeding  $\pm 2$  calibrator deviation ranges found for the controls.

**Bilirubin**

Sample #	Control Mean (pg/mL)	Control 2 SD range	Bilirubin Mean (pg/mL)	Interference ANOVA
Low	19.9	18.5-21.3	19.5	0.54
Mid	31.7	27.7-35.7	29.9	0.22
High	79.7	63.5-95.9	83.2	0.45

**Cholesterol**

Sample #	Control Mean (pg/mL)	Control 2 SD range	Cholesterol Mean (pg/mL)	Interference ANOVA
Low	19.9	18.5-21.3	19.4	0.64
Mid	31.7	27.7-35.7	29.7	0.15
High	79.7	63.5-95.9	75.2	0.40

**Triglyceride**

Sample #	Control Mean (pg/mL)	Control 2 SD range	Triglyceride Mean (pg/mL)	Interference ANOVA
Low	14.7	12.9-16.5	15.1	0.88
Mid	30.4	22.4-38.4	33.4	0.25
High	92.6	69.5-116.9	95.4	0.68

**Hemoglobin**

Sample #	Control Mean (pg/mL)	Control 2 SD range	Hemoglobin Mean (pg/mL)	Interference ANOVA
Low	21.5	17.5-25.5	20.8	0.53
Mid	38.6	35.4-41.8	40.0	0.43
High	112.7	101.4-124.0	120.6	0.17

**Urea**

Sample #	Control Mean (pg/mL)	Control 2 SD range	Urea Mean (pg/mL)	Interference ANOVA
Low	25.1	20.9-29.3	26.6	0.39
Mid	36.2	23.8-48.6	35.9	0.96
High	120.1	97.4-142.8	118.6	0.82

**REFER TO LAST PAGE FOR REFERENCES**

### SCHEME OF THE ASSAY

1. The calibrators, controls and unknowns are all reconstituted with 95% ethanol and tracer after the column procedure dry down step.
2. Dispense reagents according to the following scheme:

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Reconstituted Calibrators	-	-	75 µL	-
Reconstituted Controls and Unknown Samples	-	-	-	75 µL
95% Ethanol and Tracer **	75 µL	75 µL	-	-
NSB Buffer	-	300 µL	-	-
1,25-(OH) <sub>2</sub> -D Antiserum	-	-	300 µL	300 µL

NOTE: Total count and NSB tubes are created by adding 50 µL of 95% ethanol with 125 µL of tracer and pipetting 75 µL of this mixture into duplicate assay tubes.

3. Mix well; incubate for 2 hours (+/- 15 minutes) at 20-25°C.
4. Dispense 500 µL of GAR precipitating reagent into all wells, except the total count tubes.
5. Mix well; incubate for 20 minutes (+/- 5 minutes) at 20-25°C.
6. Centrifuge using 1800 x g\* for 20 minutes at 20-25°C.
7. Decant the supernatants.
8. Count each tube in a gamma counter for 1 minute.

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$$*g = (1118 \times 10^{-8}) (\text{radie i cm}) (\text{rpm})^2$$

**REFERENCES/BIBLOGRAPHIE/LITERATUR/BIBLIOGRAFÍA/BIBLIOGRAFIA/  
REFERENSER**


1. Haddad, J.G. and T.C.B. Stamp, "Circulating 25-Hydroxyvitamin D in Man," **American Journal of Medicine**, 57: 57, (1974).
2. Holick, M.F., "The Cutaneous Photosynthesis of Previtamin D<sub>3</sub>: A Unique Photo-endocrine System," **Journal of Investigative Dermatology**, 76: 51, (1981).
3. Harrison, J.E., A.J.W. Hitchman, G. Jones, C.S. Tam and J.N.M. Heersche, "Plasma Vitamin D Metabolite Levels in Phosphorus Deficient Rats During the Development of Vitamin D Deficient Rickets," **Metabolism**, 31 (11): 1121, (1982).
4. DeLuca, H.F., "The Vitamin D Systems: A View From Basic Science to the Clinic," **Clinical Biochemistry**, 14: 213, (1981).
5. Norman, A.W. and F.P. Roos, "Vitamin D Seco-Steroids: Unique Molecules with Both Hormone and Possible Membranophilic Properties," **Life Science**, 24: 759, (1979).
6. DeLuca, H.F. and H.K. Schnoes, "Metabolism and Mechanism of Action of Vitamin D," **Annual Reviews Biochemistry**, 45: 631, (1976).
7. Aarskog, D., L. Aksnes and T. Markestad, "Effect of Parathyroid Hormone on Vitamin D Metabolism in Osteoporosis," **Pediatrics**, 68 (1): 109, (1981).
8. Chesney, R.W., "Current Clinical Applications of Vitamin D Metabolite Research," **Clinical Orthopedics and Related Research**, 161: 285, (1981).
9. Sorensen, O.H., B. Lumholtz, B. Lund, I. Hjelmstrand, L. Mosekilde, F. Melsen, J.E. Bishop and A.W. Norman, "Acute Effects of Parathyroid Hormone on Vitamin D Metabolism in Patients with the Bone Loss of Aging," **Journal of Clinical Endocrinology and Metabolism**, 54 (6): 1258, (1982).
10. Hollis, B.W., "Assay of Circulating 1,25-dihydroxyvitamin D Involving a Novel Single-Cartridge Extraction and Purification Procedure," **Clinical Chemistry**, 32 (11): 2060, (1986).
11. Lissner, D., R.S. Mason and S. Posen, "Stability of Vitamin D Metabolites in Human Blood Serum and Plasma," **Clinical Chemistry**, 27 (5): 773, (1981).
12. Tonks, D.B., "A Study of the Accuracy and Precision of Clinical Chemistry Determinations in 170 Canadian Laboratories," **Clinical Chemistry**, 9: 217, (1963).
13. Rodbard, D., et. al., "Statistical Quality Control of Radioimmunoassays," **The Journal of Clinical Endocrinology Metabolism**, 28: 1412, (1968).
14. National Committee for Clinical Laboratory Standards, Internal Quality Control Testing: Principles and Definitions, Approved Guideline, **NCCLS** document C24-A (ISBN 1-56238-112-1), NCCLS, 771 East Lancaster Avenue, Villanova, PA, 19085, (1991).
15. National Committee for Clinical Laboratory Standards, "How to Define and Determine Reference Intervals in the Clinical Laboratory, 2<sup>nd</sup> ed." **NCCLS** Document C28-A2, 20(13) Approved Guideline, (2000).
16. National Committee for Clinical Laboratory Standards, "Evaluation of Precision Performance of Clinical Chemistry Devices," **NCCLS** Document EP5-A, 19(2), Approved Guideline, (1999).
17. National Committee for Clinical Laboratory Standards, "Evaluation of the Linearity of Quantitative Analytical Methods: Performance of Clinical Chemistry Devices," **NCCLS** Document EP6-P, 6(18), Proposed Guideline, (1996).
18. National Committee for Clinical Laboratory Standards, "Interference Testing in Clinical Chemistry," **NCCLS** Document EP7-P, 6(13), Proposed Guideline, (1986).

### SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
<b>IVD</b>	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
<b>LOT</b>	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Maximum temperature	Température maximale	Höchsttemperatur	Temperatura máxima	Temperatura massima
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
<b>Ab</b>	Antiserum	Antisérum	Antiserum	Antisero	Antisiero
<b>Ab PEG</b>	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reactivo precipitante	Reagente precipitante
<b>Ag</b> <sup>125</sup> I	Tracer: antigen labelled with <sup>125</sup> I	Traceur : antigène marqué à l'iode <sup>125</sup>	Tracer: <sup>125</sup> I-markiertes Antigen	Trazador: antígeno etiquetado con <sup>125</sup> I	Tracciatore: antigene etichettato con <sup>125</sup> I
<b>BUF</b>	Buffer	Tampon	Puffer	Tampón	Tampone
<b>CAL</b>	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
<b>CONTROL</b>	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
<b>REAG</b>	Reagent	Réactif	Reagenz	Reactivo	Reagente
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo
	Highly flammable	Facilement inflammable	Leichtentzündlich	Fácilment inflamable	Facilment infiammabile
	Harmful	Nocif	Mindergiftig	Nocivo	Nocivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,  
13040 Saluggia (VC) Italy

Phone: 651-439-9710  
FAX: 651-351-5669

For Customer Service in the  
U.S. and Canada Call Toll Free:  
800-328-1482

In the United Kingdom Call:  
44 118 9364200  
FAX: 44 118 9792061

13161E

27735  
1/04

PRINTED IN U.S.A.

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DiaSorin Inc.